

# SONNET 2 Audio Processor

## Product Catalogue



hearLIFE



# SONNET 2

## Audio Processor

### SONNET 2

SONNET 2 is the latest audio processor from MED-EL. Featuring advanced hearing technology, it is designed to give optimal sound quality. But not only that, SONNET 2 can also be customised to suit the personal needs of each individual user. From flexible battery options, to universal connectivity accessories, SONNET 2 is made for you.

### SONNET 2 EAS


SONNET 2 EAS contains Electric Acoustic Stimulation technology to assist users with high-frequency hearing loss. It combines electrical CI stimulation with acoustic amplification, to take advantage of a user's natural residual hearing. SONNET 2 EAS offers all the same great features as SONNET 2, plus an additional acoustic unit and earpiece.



# SONNET 2 User Kit

The SONNET 2 User Kit is the must-have product for all new users. From a fully-assembled audio processor, to waterproof covers and other accessories, it provides everything needed to get started with SONNET 2. What's more it all comes in a stylish backpack, which can be used either to store SONNET 2, or just for everyday use.

## How To Choose A User Kit



In addition to the fully-assembled audio processor, all user kits contain the following accessories:

- WaterWear: see page 27
- FineTuner Echo Kit: see page 14
- Daily Case Kit: see page 15
- Drying Kit: see page 16
- MED-EL backpack: see page 16
- Important documents and user manuals

Depending on which battery option is selected, the user kit will also contain the following items:

#### Rechargeable Batteries:

- 2 Rechargeable Batteries Standard
- 1 Rechargeable Battery Micro
- Rechargeable Battery Adapter
- Charging System
- Battery Wallet

#### Disposable Batteries:

- Battery Pack Frame
- Pack of Size 675 Zinc-Air Batteries

#### Both Types of Batteries:

- All of the above

A Spare Parts Kit should also be ordered alongside the user kit.

The Spare Parts Kits contain audio processor parts in different colours, so your patients can customise their SONNET 2. Select the right Spare Parts Kit by using the diagram on page 4.

## Step 01 Audio Processor

Choose the type of audio processor

01



**SONNET 2**  
CI Audio Processor



**SONNET 2 EAS**  
EAS Audio Processor

## Step 02 Implant

Choose the type of implant\*

02

**SYNCHRONY Implants**

**Other Implants**

**SYNCHRONY Implants**

**Other Implants**

## Step 03 Battery & Colour

Choose battery and colour options

03

Disposable Battery

Rechargeable Battery

Both

Disposable Battery

Rechargeable Battery

Both

Disposable Battery

Rechargeable Battery

Both

Disposable Battery

Rechargeable Battery

Both

	SYNCHRONY Implants			Other Implants			SYNCHRONY Implants			Other Implants		
	Disposable Battery	Rechargeable Battery	Both	Disposable Battery	Rechargeable Battery	Both	Disposable Battery	Rechargeable Battery	Both	Disposable Battery	Rechargeable Battery	Both
Anthracite	37581	37582	37583	37599	37600	37601	37725	37726	37727	37743	37744	37745
Beige	37593	37594	37595	37611	37612	37613	37737	37738	37739	37755	37756	37757
Black	37587	37588	37589	37605	37606	37607	37731	37732	37733	37749	37750	37751
Ebony	37584	37585	37586	37602	37603	37604	37728	37729	37730	37746	37747	37748
Nordic Grey	37596	37597	37598	37614	37615	37616	37740	37741	37742	37758	37759	37760
White	37590	37591	37592	37608	37609	37610	37734	37735	37736	37752	37753	37754

\* The SYNCHRONY series of implants includes any MED-EL implant with SYNCHRONY in the name. Implants not in the SYNCHRONY series include CONCERTO and SONATA. For a full list of implants, please visit [go.medel.com/faq-implants](http://go.medel.com/faq-implants)

# SONNET 2 Spare Parts Kits

The Spare Parts Kit contains certain SONNET 2 parts in a wide range of colours. This allows SONNET 2 to be customised to suit the style of the user.

## Choose A Spare Parts Kit

	<b>Disposable Battery Option</b> 2 Microphone Covers 1 DL-Coil Cover Low Profile 1 Battery Pack Cover Standard	<b>Rechargeable Battery / Both Batteries Option</b> 2 Microphone Covers 1 DL-Coil Cover Low Profile 1 Battery Pack Cover Standard 1 Battery Pack Cover Micro
Anthracite	37408	37420
Baby Blue	37414	37426
Baby Pink	37419	37431
Beige	37410	37422
Black	37412	37424
Bordeaux Red	37415	37427
Ebony	37409	37421
Green	37418	37430
Nordic Grey	37413	37425
Orange	37417	37429
Pacific Blue	37416	37428
White	37411	37423









# SONNET 2 Main Components

Because of its flexible modular design, SONNET 2 is easy to customise with different parts. All parts listed below can be ordered separately as required.



## SONNET 2 Control Unit







For CI and ABI—Includes Earhook (Regular), Microphone Cover, and Earhook Pin. The audio processor will need to be refitted if this component is replaced.

	Anthracite	34917		Ebony	34925
	Beige	34949		Nordic Grey	34957
	Black	34933		White	34941



## SONNET 2 EAS Control Unit
















For EAS—Includes EAS Earhook (Regular), Microphone Cover, and Earhook Pin. The audio processor will need to be refitted if this component is replaced.


	Anthracite	34921		Ebony	34929
	Beige	34953		Nordic Grey	34961
	Black	34937		White	34945



## Microphone Cover

Exchangeable covers available in 15 glossy colours. The cover should be replaced every 3 months.

	Anthracite	30719		Ebony	30858
	Baby Blue	30870		Green	30868
	Baby Pink	30869		Nordic Grey	30862
	Beige	30861		Orange	30867
	Black	30859		Pacific Blue	30866
	Bordeaux Red	30865		Sienna Brown	30863
	Creme	30864		White	30860
	Cool Grey	30871			



## CI Earhooks

CI Earhook (Regular)	30118
CI Earhook (Large)	30422
Earhook Pin	30098

## EAS Earhooks

EAS Earhook (Regular)	30720
EAS Earhook (Large)	30653
Earhook Pin	30098



# DL-Coil

SONNET 2's DL-Coil has several features, designed to give reliable hearing. The DL-Coil features a link-check indicator light, locking cable connection, durable coil cables, and adjustable coil magnets.



## DL-Coil Base





Can be used with any exchangeable DL-Coil cover and DL-Coil magnet.

	Nordic Grey	31305
---	-------------	-------



## DL-Coil Cable






Connects the DL-Coil to SONNET 2. Available in 3 lengths and 4 colours.

	6.5 cm (2.6 in.)	9 cm (3.5 in.)	28 cm (11 in.)	
	Beige	32544	32545	31582
	Black	32540	32541	31580
	Nordic Grey	32546	32547	31583
	White	32542	32543	31581



## DL-Coil Magnet SYNCHRONY

For SYNCHRONY implants only\*






	DL-Coil Magnet 15—Soft	31367
	DL-Coil Magnet 25—Medium	31368
	DL-Coil Magnet 35—Strong	31369
	DL-Coil Magnet 45—Super Strong**	31370
	DL-Coil Magnet 55—Max Strength**	31371

\*\* Use with High Profile Cover only



## DL-Coil Magnet













For all other implants\*

	DL-Coil Magnet 1—Soft	31362
	DL-Coil Magnet 2—Medium	31363
	DL-Coil Magnet 3—Strong	31364
	DL-Coil Magnet 4—Super Strong**	31365
	DL-Coil Magnet 5—Max Strength**	31366



## Low Profile







For use with Soft, Medium, or Strong DL-Coil Magnets.

	Anthracite	31159		Ebony	31160
	Baby Blue	32979		Green	32982
	Baby Pink	32980		Nordic Grey	31164
	Beige	31163		Orange	32981
	Black	31161		Pacific Blue	32978
	Bordeaux Red	32977		White	31162



## High Profile

For use with Super Strong or Max Strength DL-Coil magnets.

	Anthracite	31165		Ebony	31166
	Beige	31169		Nordic Grey	31170
	Black	31167		White	31168

\* The SYNCHRONY series of implants includes any MED-EL implant with SYNCHRONY in the name. Implants not in the SYNCHRONY series include CONCERTO and SONATA. For a full list of implants, please visit [go.medel.com/faq-implants](http://go.medel.com/faq-implants)

# Design Covers for DL-Coil

SONNET 2 is designed to blend in with your natural look. The Design Covers for DL-Coil allow the coil to be matched to specific hair tones, for those who like to keep their audio processor discreet. Alternatively, the art, animal, and pattern inspired covers enable recipients to make a statement with their SONNET 2.



Captain's Choice

Low Profile  
High Profile

35213  
35237



Superhero

Low Profile  
High Profile

35211  
35236



Green Monster

Low Profile  
High Profile

35115  
35233



Princess

Low Profile  
High Profile

35210  
35235



Kitten

Low Profile  
High Profile

35116  
35234



Stained Glass

Low Profile  
High Profile

35217  
35241



Melody

Low Profile  
High Profile

35216  
35240



Coral

Low Profile  
High Profile

35218  
35242



Leopard Skin

Low Profile  
High Profile

35214  
35238



Zebra Skin

Low Profile  
High Profile

35215  
35239



Onyx

Low Profile  
High Profile

35220  
35244



Oriental Flower

Low Profile  
High Profile

35223  
35247



Tartan

Low Profile  
High Profile

35221  
35245



Watercolour

Low Profile  
High Profile

35222  
35246



Wood

Low Profile  
High Profile

35219  
35243



Silver Grey

Low Profile  
High Profile

35224  
35248



Platinum Blonde

Low Profile  
High Profile

35225  
35249



Blonde

Low Profile  
High Profile

35226  
35250



Straw Blonde

Low Profile  
High Profile

35227  
35251



Light Brown

Low Profile  
High Profile

35228  
35252



Golden Brown

Low Profile  
High Profile

35229  
35253



Auburn

Low Profile  
High Profile

35230  
35254



Dark Brown

Low Profile  
High Profile

35231  
35255



Jet Black

Low Profile  
High Profile

35232  
35256



## SONNET 2 Battery Options

Choose from rechargeable batteries, disposable batteries, or both! SONNET 2's flexible battery options allow users to power their audio processor how it best suits them.

In order to use both type of batteries, users must have both the Rechargeable Battery Adapter (for rechargeable batteries) and the Battery Pack Frame (for disposable batteries).

### Rechargeable Batteries



#### Rechargeable Battery Standard

Lithium-ion battery provides up to 10 hours of use per charge.

Pack with 1 Battery  
 Pack with 2 Batteries  
 Pack with 3 Batteries

31267  
 33634  
 33632



#### Rechargeable Battery Micro

Compact lithium-ion battery provides up to 7 hours of use per charge.

Pack with 1 Battery  
 Pack with 2 Batteries  
 Pack with 3 Batteries


31266  
 33633  
 34143



#### Rechargeable Battery Adapter

Connects a rechargeable battery pack to the SONNET 2 Control Unit.

32218



**Battery Wallet**

Battery carrying wallet with three pockets for both charged and empty batteries.

33719



**Charging System**

33643

USB-powered charger can charge up to 4 rechargeable SONNET 2 batteries at once.

- Charging Unit
- USB Cable
- Soft Carrying Pouch

- 32244
- 37022
- 33889



**USB Power Supply**

37828

USB power supply for Charging System. Includes adapter for US, UK, EU, and AU.

# Disposable Batteries



Battery Pack Frame

For use with 2 size 675 batteries.

30467



Disposable Batteries

Pack of 6 mercury-free, high-power zinc-air size 675 batteries.

09564

# Battery Pack Covers



Battery Pack Cover Standard

For use with the Rechargeable Battery Standard, and the Battery Pack Frame with disposable batteries. It can be used with WaterWear.

Anthracite	30881	Ebony	30788
Baby Blue	30879	Green	30877
Baby Pink	30878	Nordic Grey	30792
Beige	30791	Orange	30876
Black	30789	Pacific Blue	30875
Bordeaux Red	30874	Sienna Brown	30872
Creme	30873	White	30790
Cool Grey	30880		



Battery Pack Cover Micro

For use with the Rechargeable Battery Micro. It cannot be used with WaterWear.

Anthracite	32219	Ebony	32220
Baby Blue	32232	Green	32230
Baby Pink	32231	Nordic Grey	32224
Beige	32223	Orange	32229
Black	32221	Pacific Blue	32228
Bordeaux Red	32227	Sienna Brown	32225
Creme	32226	White	32222
Cool Grey	32233		




# Mini Battery Pack for SONNET 2



The Mini Battery Pack powers SONNET 2 for up to 37 hours on a single AAA battery. It can also run for up to 16 hours on a single rechargeable DaCapo PowerPack.



## Mini Battery Pack Kit, 50 cm




1 Mini Battery Pack Cable (50 cm)	
1 Mini Battery Pack (Black)	05880
1 Mini Battery Pack Adapter	32290
1 AAA Battery	08083
1 Mini Battery Pack Accessories	32579



	Anthracite Kit	33372
	Cable	32291
	Beige Kit	33377
	Cable	32295
	Black Kit	33375
	Cable	32293

	Ebony Kit	33374
	Cable	32292
	Nordic Grey Kit	33378
	Cable	32296
	White Kit	33376
	Cable	32294

## Mini Battery Pack Kit, 75 cm

1 Mini Battery Pack Cable (75 cm)	
1 Mini Battery Pack (Black)	05880
1 Mini Battery Pack Adapter	32290
1 AAA Battery	08083
1 Mini Battery Pack Accessories	32579

	Anthracite Kit	33379
	Cable	32297
	Beige Kit	33383
	Cable	32301
	Black Kit	33381
	Cable	32299

	Ebony Kit	33380
	Cable	32298
	Nordic Grey Kit	33384
	Cable	32302
	White Kit	33382
	Cable	32300

# FineTuner Echo

The FineTuner Echo makes it simple to change the settings on SONNET 2. With an easy-to-read display, plus integrated microphone check function, the FineTuner Echo is a remote control designed with the user in mind.

## FineTuner Echo Kit

Contains one FineTuner Echo including battery, plus the FineTuner Accessories Kit.

37572



## FineTuner Echo

CR2032 Battery  
Battery Lid  
Form Ring

37342

05745  
05904  
05907

## FineTuner Accessories Kit

1 Lanyard  
1 Display Foil  
1 Microfibre Cloth

37343

08248  
36232  
35915





## Daily Case Kits

The Daily Case is the perfect accessory for carrying and storing SONNET 2. Each kit also contains a spare earhook, plus tools for cleaning your audio processor.

### Daily Case Kit

1 Daily Case M	34706
1 Brush	05097
1 Audio Processor Combi Tool	07486
1 Earhook for CI or	30422
1 Earhook for EAS	30653
1 Earhook Pin	30098
1 Pack Zinc-Air Batteries	09564

	For CI	For EAS
With Disposable Batteries	37568	37569
Without Disposable Batteries	37570	37571





# SONNET 2 Accessories

These convenient accessories are ideal for everyday life with your SONNET 2.



## All-Purpose Backpack

High-quality backpack for carrying accessories and other items. 35280



## Dry-Space UV

Daily drying helps keep SONNET 2 in optimal working condition, and UV light sanitises the audio processors. Each drying kit has space for two SONNET 2 processors. Powered by included AC adapter. 36722

## Audio Processor Travel Case

Convenient zippered carrying case for audio processor and accessories. 33722

Daily Case M 34706

Daily Case S 34705



# AudioLink

AudioLink is a universal connectivity device that can be used to connect SONNET 2 to phones, tablets, TVs, and many other media devices.

<b>AudioLink Kit*</b>	<b>37218</b>
1 AudioLink Main Unit	37028
1 Docking Station	37034
1 Opening Tool	36795
1 Cable Kit	37824
1 Accessories Kit	37825
*including Battery Unit	



<b>AudioLink Spare Battery Kit</b>	<b>38058</b>
1 Battery Unit	37857
1 Opening Tool	36795



<b>AudioLink Accessories Kit</b>	<b>37825</b>
1 Clip	34307
1 Lanyard	08248



<b>AudioLink Cable Kit</b>	<b>37824</b>
1 Micro USB Cable	37022
1 RCA Cable	35512
1 Audio Cable	36228
1 Audio Cable Adapter	36792





# Neckloops

Bluetooth neckloops are a simple and convenient way to stream audio wirelessly to SONNET 2. You can connect to your phone through Bluetooth or plug into your favourite hand-held devices.



Artone 3 MAX Neckloop

Integrated Bluetooth and dual microphones enable hands-free phone calls and music streaming.

31776



Geemarc iLoop+ Neckloop

Connects to any standard audio headphone jack with a 3.5 mm plug.

08320









# Direct Audio Input Cables

Connect SONNET 2 to audio devices, like phones, music players, and tablets. The 3.5 mm jack at one end of the cable plugs into any standard headphone port. The 3-pin jack at the other end connects the SONNET 2 FM Battery Pack Cover.



## FM Battery Pack Cover

Enables connection to 3-pin FM systems or direct audio input cables.

	Anthracite	30492		Ebony	30833
	Beige	30836		Nordic Grey	30837
	Black	30834		White	30835



## Unilateral Mixed Source

Provides a 50/50 mix of sound from the audio device and the SONNET 2 microphones, depending on settings. For use with one audio processor. 04438



## Bilateral Mixed Source

Provides a 50/50 mix of sound from the audio device and the SONNET 2 microphones, depending on settings. For use with two audio processors. 04440



## Unilateral External Source

Provides a 90/10 mix of sound, with more sound from the audio device, depending on settings. For use with one audio processor. 04439



## Bilateral External Source

Provides a 90/10 mix of sound, with more sound from the audio device, depending on settings. For use with two audio processors. 04441



# Wearing Options

## ActiveWear for SONNET 2

ActiveWear is the smallest audio processor wearing option available and is the ideal choice for younger children, or active adults. It allows the battery pack to be attached to clothing, leaving just the lightweight control unit behind the ear. With ActiveWear, children have a lighter way to wear SONNET 2 all day long, whereas adults have a more secure option during sports.



### SONNET 2 ActiveWear Kit, 27cm

- 1 ActiveWear Cable (27 cm)
- 1 Design Cover - Melody 35216
- 1 Design Cover - Onyx 35220
- 1 Attachment Clip 35040
- 1 Huggy 02547
- 1 Activity Hook 05751

Anthracite Kit	35124
Anthracite Cable	34594
Beige Kit	35128
Beige Cable	34598
Black Kit	35126
Black Cable	34596

Ebony Kit	35125
Ebony Cable	34595
Nordic Grey Kit	35129
Nordic Grey Cable	34599
White Kit	35127
White Cable	34597

### SONNET 2 ActiveWear Kit, 50cm

- 1 ActiveWear Cable (50 cm)
- 1 Design Cover - Melody 35216
- 1 Design Cover - Onyx 35220
- 1 Attachment Clip 35040
- 1 Huggy 02547
- 1 Activity Hook 05751

Anthracite Kit	35130
Anthracite Cable	34600
Beige Kit	35134
Beige Cable	34604
Black Kit	35132
Black Cable	34602

Ebony Kit	35131
Ebony Cable	34601
Nordic Grey Kit	35135
Nordic Grey Cable	34605
White Kit	35133
White Cable	34603





## BabyWear for SONNET 2

BabyWear is a comfortable wearing option for babies and toddlers who aren't yet ready to wear a behind-the-ear processor. With BabyWear, SONNET 2 is worn completely off the ear and is instead attached to a baby's clothing. Alternatively, BabyWear can be used by children and adults when wearing a helmet, or in any other situation when the audio processor can't be worn on the ear.



### SONNET 2 BabyWear Kit for DL-Coil

1 DL-Coil Cable (28 cm)	
1 Design Cover - Green Monster	35115
1 Design Cover - Kitten	35116
1 Earhook Cover Set	35499
1 Attachment Clip	35040
1 Huggy	02547

Beige Kit	35122	Nordic Grey Kit	35123
Black Kit	35120	White Kit	35121

# Additional Fixation

These easy-to-use accessories help keep SONNET 2 right where it should be. They are ideal for young children or for additional security during activities.



BTE Headband

	S (43-47 cm)	M (48-51cm)
<span style="color: pink;">■</span> Pink	33847	33848
<span style="color: blue;">■</span> Blue	33839	33846



Mellie Clip

Lanyard that connects to an audio processor and can be clipped onto clothing.

02295



Activity Hook

Slides onto SONNET 2 to provide additional support on the bottom of the ear.

05751



Huggy

Soft fixation that goes over the full ear to hold the audio processor in place or attach it to clothing.

02547



GripWear

Adhesive pad which is applied to the bottom of the DL-Coil for better fixation. Includes 12 pads per package.

33686





## WaterWear for SONNET 2

WaterWear is a fully waterproof accessory cover for SONNET 2 that can be used in fresh-, salt-, or chlorinated water. Each pack contains 3 reusable covers, plus adhesive strips, for multiple trips to the pool or beach.

Only the standard SONNET 2 Battery Pack Cover with rechargeable, alkaline, or silver-oxide batteries can be used with WaterWear.



### WaterWear for SONNET 2

---

Includes 3 reusable covers, plus adhesive strips.

38232

### Batteries for WaterWear

---

Alkaline Batteries LR44

32328



## Get In Touch

### Where to Order

To order parts and accessories for SONNET 2, contact your local MED-EL representative.

Let us know if you have any questions—we look forward to hearing from you!



[medel.com](http://medel.com)



## MED-EL Offices Worldwide

### AMERICAS

**Argentina**  
medel@medel.com.ar

**Brazil**  
office@br.medel.com

**Canada**  
officecanada@medel.com

**Colombia**  
office-colombia@medel.com

**Mexico**  
office-mexico@medel.com

**United States**  
implants@medelus.com

### ASIA PACIFIC

**Australia**  
office@medel.com.au

**China**  
office@medel.net.cn

**Hong Kong**  
office@hk.medel.com

**India**  
implants@medel.in

**Indonesia**  
office@id.medel.com

**Japan**  
office-japan@medel.com

**Malaysia**  
office@my.medel.com

**Philippines**  
office@ph.medel.com

**Singapore**  
office@sg.medel.com

**South Korea**  
office@kr.medel.com

**Thailand**  
office@th.medel.com

**Vietnam**  
office@vn.medel.com

### EMEA

**Austria**  
office@at.medel.com

**Belgium**  
office@be.medel.com

**Denmark**  
nordic@medel.com

**Finland**  
office@fi.medel.com

**France**  
office@fr.medel.com

**Germany**  
office@medel.de

**Italy**  
ufficio.italia@medel.com

**Netherlands**  
office@nl.medel.com

**Norway**  
norge@medel.com

**Portugal**  
office@pt.medel.com

**South Africa**  
customerserviceZA@medel.com

**Spain**  
office@es.medel.com

**Sweden**  
nordic@medel.com

**Switzerland**  
office-ch@medel.com

**Turkey**  
medelturkey@medel.com

**United Arab Emirates**  
office@ae.medel.com

**United Kingdom**  
customerservice@medel.co.uk