

OPUS 2

Patient Kit Product Catalogue





OPUS 2 Basekit

07584

OPUS 2 Base Kit Anthracite* with D Coil, consisting of:

	04310 OPUS 2 Control Unit Anthracite*		04313 OPUS 2 Battery Pack Anthracite*
	06131 D Coil Anthracite*		06121 2 x Coil Cable 7.5 cm for D Coil
	04273 Connecting Piece Anthracite*		04449 2 x Connecting Piece for Safety Lock Anthracite*
	04463 2 x Safety Lock (in combination with Connecting Piece)		04048 Speech Processor Test Device (grey)
	04278 OPUS 2 Earhook (Spare)		04454 OPUS 2 Earhook Pin (Spare)
	03132 OPUS 2 FineTuner		02910 2 x 6 Zinc Air Batteries (Type 675)
	03142 CR2025 Battery (included in FineTuner)		05097 Cleaning Brush

OPUS 2 Basekit (cont.)



04653
Drying Kit



05862
Electrical Drying Kit



03171
Drying Kit Tablets



04637
Key Ring (for storage of
zinc air batteries)



05093
Screwdriver



05094
Pin Removal Tool

04724
OPUS 2 Patient Kit
Packaging List

04725
OPUS 2 Registration Card
(with envelope)

05266
OPUS 2 Support Kit

OPUS 2 User Manual
OPUS 2 Patient Kit Case

OPUS 2 Coloured Battery Pack Cover Extension Kit, consisting of:



04289
OPUS 2 Battery Pack
Cover Pacific Blue

04290
OPUS 2 Battery Pack
Cover Bordeaux Red

05859
OPUS 2 Battery Pack
Cover Baby Blue

05858
OPUS 2 Battery Pack
Cover Baby Pink

06655
OPUS 2 Battery Pack
Cover Orange

06656
OPUS 2 Battery Pack
Cover Green

*different colours available; please refer to page 11

OPUS 2XS Extension Kit

09242

OPUS 2XS Extension Kit anthracite*, consisting of:



08091

XS Battery Pack
anthracite*



05751

Activity Hook



04675

Battery Pack Pin
Protection Cover

XS Battery Pack
User Manual

*different colours available; please refer to page 11

OPUS 2 Children's Extension Kit

05000

OPUS 2 Children Extension Kit Anthracite*, consisting of:



02421

Children's Battery Pack Anthracite* (with 27 cm cable)



01890

Battery Pack Lid Anthracite*



04461

Children's/Mini Battery Pack Connecting Piece



04462

2 x Children's/Remote Battery Pack Connecting Piece for Safety Lock



04463

2 x Safety Lock (in combination with Connecting Piece)



03108

Fixation Clip (for fixation of Children's Battery Pack to clothes)



02558

Fixation Bar (for fixation of Children's Battery Pack to clothes)



02547

Huggy



02295

Dino Clip



04675

Battery Pack Pin Protection Cover

*different colours available; please refer to page 11

OPUS 2 BabyBTE/ActiveWear Extension Kit

07599

OPUS 2 BabyBTE/ActiveWear Extension Kit Anthracite* for D Coil, consisting of:

	01862 Straight Battery Pack Anthracite*		06122 Coil Cable 28 cm for D Coil
	02468 BabyBTE/ActiveWear Connecting Piece		04465 2 x BabyBTE/ActiveWear Connecting Piece for Safety Lock
	04463 2 x Safety Lock (in combination with Connecting Piece)		03108 Fixation Clip (for fixation of BabyBTE/ActiveWear to clothes)
	02558 Fixation Bar (for fixation of BabyBTE/ActiveWear to clothes)		01890 Battery Pack Lid Anthracite*
	04270 2 x Microphone Cover		04675 Battery Pack Pin Protection Cover

*different colours available; please refer to page 11

Mini Battery Pack Kit

08264

Mini Battery Pack Kit (28 cm), consisting of:

08265

Mini Battery Pack Kit (100 cm), consisting of:



05880

Mini Battery Pack
Anthracite



08083

AAA Battery,
non-rechargeable



04803

Mini Battery Pack Cable,
28cm (included in 08264)



04804

Mini Battery Pack Cable,
100cm (included in
08265)



08246

Mini Battery Pack Pouch



08247

Mini Battery Pack Clip



08248

Lanyard



04461

Children's/Mini Battery
Pack Connecting Piece



04463

2 × Safety Lock (in
combination with
Connecting Piece)



04462

2 × Children's/Remote
Battery Pack Connecting
Piece for Safety Lock

Mini Battery Pack
User Manual

OPUS 2 FM Extension Kit

04999

OPUS 2 FM Extension Kit Anthracite*, consisting of:



04455

OPUS 2 FM Battery Pack Cover Anthracite*



04438

Adapter Cable (3.5 mm, unilateral, MIX mode – for connection to battery powered audio devices such as portable CD players, MP3 Players)

*different colours available; please refer to page 11

OPUS 2 DaCapo System Extension Kit

05378

OPUS 2 DaCapo System Extension Kit Anthracite*, consisting of:



04563

DaCapo Housing (Battery Pack Frame and Cover) Anthracite*



04392

DaCapo Charger



03602

3 x DaCapo PowerPack (rechargeable battery)



05198

DaCapo Power Pack Inlay (to be inserted in OPUS 2 Support Kit)



04675

Begriff Battery Pack Pin Protection Cover













DaCapo User Manual

*different colours available; please refer to page 11

Accessories

	<p>04439 Adapter Cable (3.5 mm, unilateral, EXT mode – for connection to battery powered audio devices such as portable CD players, MP3 Players)</p>		<p>04440 Bilateral Adapter Cable, MIX mode (with yellow plug indicating MIX mode)</p>
	<p>04441 Bilateral Adapter Cable, EXT mode (with red plug indicating EXT mode)</p>		<p>05438 DaCapo Car Adapter – for connection of DaCapo Charger to car cigarette lighter</p>
	<p>04271 Microphone Cover „left“ (L for left side engraved)</p>		<p>04279 OPUS 2 Earhook „left“ (L for left side engraved)</p>
	<p>07789 D Coil MET (Magnet Exchange Tool)</p>		<p>05751 Activity Hook</p>
	<p>05735 Children's Connecting Piece for Huggy</p>		<p>06120 Coil Cable 12 cm for D Coil</p>

OPUS 2: Available Colours

 Anthracite	 Baby Pink
 Nordic Grey	 Bordeaux Red
White	 Orange
 Creme	 Green
 Beige	 Sienna Brown
 Baby Blue	 Ebony
 Pacific Blue	

Colours shown may vary slightly from actual product colours.

OPUS 2: Base Kit and Extension Kit Overview

07584	OPUS 2 D Coil Base Kit Anthracite
07585	OPUS 2 D Coil Base Kit Beige
07586	OPUS 2 D Coil Base Kit Bordeaux Red
07587	OPUS 2 D Coil Base Kit Pacific Blue
07588	OPUS 2 D Coil Base Kit Nordic Grey
07589	OPUS 2 D Coil Base Kit Sienna Brown
07590	OPUS 2 D Coil Base Kit White
07591	OPUS 2 D Coil Base Kit Ebony
07592	OPUS 2 D Coil Base Kit Creme
07593	OPUS 2 D Coil Base Kit Orange
07594	OPUS 2 D Coil Base Kit Green
07595	OPUS 2 D Coil Base Kit Baby Blue
07596	OPUS 2 D Coil Base Kit Baby Pink
04996	OPUS 2 Colour Battery Pack Extension Kit (included in OPUS 2 Base Kit)
09242	OPUS 2XS Extension Kit Anthracite
09243	OPUS 2XS Extension Kit Beige
09244	OPUS 2XS Extension Kit Nordic Grey
09245	OPUS 2XS Extension Kit Pacific Blue
09246	OPUS 2XS Extension Kit Bordeaux Red
09247	OPUS 2XS Extension Kit Sienna Brown
09248	OPUS 2XS Extension Kit White
09249	OPUS 2XS Extension Kit Ebony
09250	OPUS 2XS Extension Kit Creme
09251	OPUS 2XS Extension Kit Orange
09252	OPUS 2XS Extension Kit Green
09253	OPUS 2XS Extension Kit Baby Blue
09254	OPUS 2XS Extension Kit Baby Pink
05000	OPUS 2 Children's Extension Kit Anthracite
05004	OPUS 2 Children's Extension Kit Beige
05008	OPUS 2 Children's Extension Kit Bordeaux Red
05012	OPUS 2 Children's Extension Kit Pacific Blue
05016	OPUS 2 Children's Extension Kit Nordic Grey
05020	OPUS 2 Children's Extension Kit Sienna Brown
07444	OPUS 2 Children's Extension Kit White

07445	OPUS 2 Children's Extension Kit Ebony
07446	OPUS 2 Children's Extension Kit Creme
07447	OPUS 2 Children's Extension Kit Orange
07448	OPUS 2 Children's Extension Kit Green
07449	OPUS 2 Children's Extension Kit Baby Blue
07450	OPUS 2 Children's Extension Kit Baby Pink
07599	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Anthracite
07600	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Beige
07601	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Bordeaux Red
07602	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Pacific Blue
07603	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Nordic Grey
07604	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Sienna Brown
07605	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit White
07606	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Ebony
07607	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Creme
07608	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Orange
07609	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Green
07610	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Baby Blue
07611	OPUS 2 D Coil BabyBTE/Active Wear Extension Kit Baby Pink
04999	OPUS 2 FM Extension Kit Anthracite
05003	OPUS 2 FM Extension Kit Beige
05007	OPUS 2 FM Extension Kit Bordeaux Red
05011	OPUS 2 FM Extension Kit Pacific Blue
05015	OPUS 2 FM Extension Kit Nordic Grey
05019	OPUS 2 FM Extension Kit Sienna Brown
07437	OPUS 2 FM Extension Kit White
07438	OPUS 2 FM Extension Kit Ebony
07439	OPUS 2 FM Extension Kit Creme
07440	OPUS 2 FM Extension Kit Orange
07441	OPUS 2 FM Extension Kit Green
07442	OPUS 2 FM Extension Kit Baby Blue
07443	OPUS 2 FM Extension Kit Baby Pink
05378	OPUS 2 DaCapo Extension Kit Anthracite
05379	OPUS 2 DaCapo Extension Kit Beige
05380	OPUS 2 DaCapo Extension Kit Bordeaux Red
05381	OPUS 2 DaCapo Extension Kit Pacific Blue
05382	OPUS 2 DaCapo Extension Kit Nordic Grey

05383	OPUS 2 DaCapo Extension Kit Sienna Brown
07458	OPUS 2 DaCapo Extension Kit White
07459	OPUS 2 DaCapo Extension Kit Ebony
07460	OPUS 2 DaCapo Extension Kit Creme
07461	OPUS 2 DaCapo Extension Kit Orange
07462	OPUS 2 DaCapo Extension Kit Green
07463	OPUS 2 DaCapo Extension Kit Baby Blue
07464	OPUS 2 DaCapo Extension Kit Baby Pink
08264	Mini Battery Pack Kit (28cm)
08265	Mini Battery Pack Kit (100cm)



MED-EL Headquarters
Fürstenweg 77a
6020 Innsbruck, Austria
office@medel.com

medel.com

