

RONDO[®] 2 Audio Processor

Product Catalogue



RONDO 2

Audio Processor

RONDO 2 is the perfect choice for recipients who want an incredibly easy-to-use and discreet wearing option. RONDO 2 is the first and only CI audio processor with an integrated wirelessly rechargeable battery, which eliminates the hassles of changing batteries.

Just place it on the charging pad and RONDO 2 starts charging. Four hours of charging provides a full day of hearing (up to 18 hours). The Mini Battery Pack also provides a back-up option to directly power RONDO 2 with a AAA battery.



RONDO 2 Core Kits



RONDO 2 Core Kit for SYNCHRONY

The RONDO 2 Core Kit includes everything necessary to get started with RONDO 2. It does not include Design Covers. Please note, packages are subject to change and content may vary. Additional products, such as Design Covers, WaterWear, or the Mini Battery Pack Kit, need to be ordered separately.

36370

- | | |
|--------------------------------|----------------------------------|
| 1 RONDO 2 Processor Unit | 1 GripWear |
| 1 RONDO 2 Magnet 3S | 1 Daily Case M |
| 1 Charging Pad with USB Cable | 1 USB Power Supply |
| 1 FineTuner | Supporting Manuals and Documents |
| 1 Speech Processor Test Device | |
| 1 Attachment Clip Set | |

RONDO 2 Core Kit Colour Mix for SYNCHRONY

The RONDO 2 Core Kit Colour Mix includes everything necessary to get started with RONDO 2. It also includes the Design Covers Colour Mix. Please note, packages are subject to change and content may vary. Additional products, such as WaterWear or the Mini Battery Pack Kit, need to be ordered separately.

35352

- | | |
|-----------------------------------|----------------------------------|
| 1 RONDO 2 Processor Unit | 1 Attachment Clip Set |
| 1 RONDO 2 Magnet 3S | 1 GripWear |
| 1 Charging Pad with USB Cable | 1 Daily Case M |
| 1 RONDO 2 Design Cover Colour Mix | 1 USB Power Supply |
| 1 FineTuner | Supporting Manuals and Documents |
| 1 Speech Processor Test Device | |

RONDO 2 Core Kit for MAESTRO

The RONDO 2 Core Kit includes everything necessary to get started with RONDO 2. It does not include Design Covers. Please note, packages are subject to change and content may vary. Additional products, such as Design Covers, WaterWear, or the Mini Battery Pack Kit, need to be ordered separately.

36369

- | | |
|--------------------------------|----------------------------------|
| 1 RONDO 2 Processor Unit | 1 GripWear |
| 1 RONDO 2 Magnet 3 | 1 Daily Case M |
| 1 Charging Pad with USB Cable | 1 USB Power Supply |
| 1 FineTuner | Supporting Manuals and Documents |
| 1 Speech Processor Test Device | |
| 1 Attachment Clip Set | |

RONDO 2 Core Kit Colour Mix for MAESTRO

The RONDO 2 Core Kit Colour Mix includes everything necessary to get started with RONDO 2. It also includes the Design Covers Colour Mix. Please note, packages are subject to change and content may vary. Additional products, such as WaterWear or the Mini Battery Pack Kit, need to be ordered separately.

35351

- | | |
|-----------------------------------|----------------------------------|
| 1 RONDO 2 Processor Unit | 1 Attachment Clip Set |
| 1 RONDO 2 Magnet 3 | 1 GripWear |
| 1 Charging Pad with USB Cable | 1 Daily Case M |
| 1 RONDO 2 Design Cover Colour Mix | 1 USB Power Supply |
| 1 FineTuner | Supporting Manuals and Documents |
| 1 Speech Processor Test Device | |

RONDO 2 Prime Kits



RONDO 2 Prime Kit for SYNCHRONY

The RONDO 2 Prime Kit is a premium kit that includes everything from the Core Kit. It does not include Design Covers. Additionally, the Prime Kit includes a high-quality carrying case, Bluetooth neckloop, and Mini Battery Pack for RONDO 2. Please note, packages are subject to change and content may vary.

36368

- | | |
|--------------------------------|--------------------------------------|
| 1 RONDO 2 Processor Unit | 1 Daily Case M |
| 1 RONDO 2 Magnet 3S | 1 USB Power Supply |
| 1 Charging Pad with USB Cable | 1 Mini Battery Pack with Accessories |
| 1 FineTuner with Accessories | 1 Artone 3 MAX Neckloop |
| 1 Speech Processor Test Device | 1 All-Purpose Backpack |
| 1 Attachment Clip Set | Supporting Manuals and Documents |
| 1 GripWear | |

RONDO 2 Prime Kit Colour Mix for SYNCHRONY

The RONDO 2 Prime Kit Colour Mix is a premium kit that includes everything from the Core Kit. Additionally, the Prime Kit Colour Mix includes a high-quality carrying case, Bluetooth neckloop, Mini Battery Pack, and the Design Cover Colour Mix. Please note, packages are subject to change and content may vary.

35321

- | | |
|-----------------------------------|--------------------------------------|
| 1 RONDO 2 Processor Unit | 1 GripWear |
| 1 RONDO 2 Magnet 3S | 1 Daily Case M |
| 1 Charging Pad with USB Cable | 1 USB Power Supply |
| 1 RONDO 2 Design Cover Colour Mix | 1 Mini Battery Pack with Accessories |
| 1 FineTuner with Accessories | 1 Artone 3 MAX Neckloop |
| 1 Speech Processor Test Device | 1 All-Purpose Backpack |
| 1 Attachment Clip Set | Supporting Manuals and Documents |

RONDO 2 Prime Kit for MAESTRO

The RONDO 2 Prime Kit is a premium kit that includes everything from the Core Kit. It does not include Design Covers. Additionally, the Prime Kit includes a high-quality carrying case. Please note, packages are subject to change and content may vary.

36367

- | | |
|--------------------------------|--------------------------------------|
| 1 RONDO 2 Processor Unit | 1 Daily Case M |
| 1 RONDO 2 Magnet 3 | 1 USB Power Supply |
| 1 Charging Pad with USB Cable | 1 Mini Battery Pack with Accessories |
| 1 FineTuner with Accessories | 1 Artone 3 MAX Neckloop |
| 1 Speech Processor Test Device | 1 All-Purpose Backpack |
| 1 Attachment Clip Set | Supporting Manuals and Documents |
| 1 GripWear | |

RONDO 2 Prime Kit Colour Mix for MAESTRO

The RONDO 2 Prime Kit Colour Mix is a premium kit that includes everything from the Core Kit. Additionally, the Prime Kit Colour Mix includes a high-quality carrying case, and the Design Cover Colour Mix for RONDO 2. Please note, packages are subject to change and content may vary.

35320

- | | |
|-----------------------------------|--------------------------------------|
| 1 RONDO 2 Processor Unit | 1 GripWear |
| 1 RONDO 2 Magnet 3 | 1 Daily Case M |
| 1 Charging Pad with USB Cable | 1 USB Power Supply |
| 1 RONDO 2 Design Cover Colour Mix | 1 Mini Battery Pack with Accessories |
| 1 FineTuner with Accessories | 1 Artone 3 MAX Neckloop |
| 1 Speech Processor Test Device | 1 All-Purpose Backpack |
| 1 Attachment Clip Set | Supporting Manuals and Documents |

RONDO 2 Main Components



RONDO 2 Processor Unit

RONDO 2 Audio Processor with integrated coil and rechargeable battery. Please note, a magnet is not included and should be ordered separately.

30459



Rubber Socket Cover

Replacement seal for RONDO 2 auxiliary input port.

33019



Magnet SYNCHRONY

For SYNCHRONY implant only.

▲	1S—Soft	35153
▲▲	2S—Medium	35154
▲▲▲	3S—Strong	35155



Magnet

For all implants except SYNCHRONY.

●	1—Soft	31253
●●	2—Medium	31254
●●●	3—Strong	31255
●●●●	4—Super Strong	31256



Magnet Exchange Tool SUP

For clinicians to exchange RONDO 2 magnets.

33422



RONDO 2 Covers

The RONDO 2 Covers come in several neutral shades, and are quick and easy to change. To maintain best possible hearing with RONDO 2, the cover should be replaced every 3 months.



Matt Cover

Single Cover

	Anthracite	35547
	Black	35546
	Cool Grey	35548
	Ebony	35550
	Nordic Grey	35549

Pack with 4 Covers

	Anthracite	36059
	Black	36058
	Cool Grey	36060
	Ebony	36062
	Nordic Grey	36061



Glossy Cover

Single Cover

	Anthracite	35553
	Black	35552
	Cool Grey	35554
	Ebony	35556
	Nordic Grey	35555

Pack with 4 Covers

	Anthracite	36064
	Black	36063
	Cool Grey	36065
	Ebony	36067
	Nordic Grey	36066

Design Covers for RONDO 2

RONDO 2 is designed to blend in with your natural look. The Design Covers allow the processor to be matched to specific hair tones, for those who like to keep their audio processor discreet. Alternatively, the art, animal, and pattern inspired covers enable recipients to show off their style with their RONDO 2.



Captain's Choice

Single Cover
Pack with 4 Covers

35659
36068



Superhero

Single Cover
Pack with 4 Covers

35658
36069



Green Monster

Single Cover
Pack with 4 Covers

35655
36070



Princess

Single Cover
Pack with 4 Covers

35657
36071



Kitten

Single Cover
Pack with 4 Covers

35656
36072



Stained Glass

Single Cover
Pack with 4 Covers

35663
36073



Melody

Single Cover
Pack with 4 Covers

35662
36074



Coral

Single Cover
Pack with 4 Covers

35664
36075



Leopard Skin

Single Cover
Pack with 4 Covers


35660
36076



Zebra Skin

Single Cover
Pack with 4 Covers

35661
36077



Onyx

Single Cover
Pack with 4 Covers

35666
36078



Oriental Flower

Single Cover
Pack with 4 Covers

35669
36079



Tartan

Single Cover
Pack with 4 Covers

35667
36080



Watercolour

Single Cover
Pack with 4 Covers

35668
36081



Wood

Single Cover
Pack with 4 Covers

35665
36082



Silver Grey

Single Cover
Pack with 4 Covers

35670
36083



Platinum Blonde

Single Cover
Pack with 4 Covers

35671
36084



Blonde

Single Cover
Pack with 4 Covers

35672
36085



Straw Blonde

Single Cover
Pack with 4 Covers

35673
36086



Light Brown

Single Cover
Pack with 4 Covers

35674
36087



Golden Brown

Single Cover
Pack with 4 Covers

35675
36088



Auburn

Single Cover
Pack with 4 Covers

35676
36089



Dark Brown

Single Cover
Pack with 4 Covers

35677
36090



Jet Black

Single Cover
Pack with 4 Covers

35678
36091

RONDO 2 Cover Kits



Black Cover Kit

Set of 4 Covers 36371

- 1 Matt Black Cover
- 1 Glossy Black Cover
- 1 Jet Black Design Cover
- 1 Watercolour Design Cover



Anthracite Cover Kit

Set of 4 Covers 36372

- 1 Matt Anthracite Cover
- 1 Glossy Anthracite Cover
- 1 Jet Black Design Cover
- 1 Watercolour Design Cover



Ebony Cover Kit

Set of 4 Covers 36375

- 1 Matt Ebony Cover
- 1 Glossy Ebony Cover
- 1 Dark Brown Design Cover
- 1 Watercolour Design Cover



Nordic Grey Cover Kit

Set of 4 Covers 36374

- 1 Matt Nordic Grey Cover
- 1 Glossy Nordic Grey Cover
- 1 Light Brown Design Cover
- 1 Watercolour Design Cover



Cool Grey Cover Kit

Set of 4 Covers 36373

- 1 Matt Cool Grey Cover
- 1 Glossy Cool Grey Cover
- 1 Silver Grey Design Cover
- 1 Watercolour Design Cover

RONDO 2 Additional Components



FineTuner Remote Control

The FineTuner remote is used to change settings such as volume, program, and connectivity settings. FineTuner is compatible with all current MED-EL audio processors. One FineTuner can be used to control both processors for bilateral users. Includes 1 CR2025 battery.

03132



Charging Pad

Wireless charging pad for RONDO 2. Powered by USB. Charging Pad not manufactured by MED-EL.

35170

FineTuner Battery

CR2025

03142

USB Power Supply

USB power supply for RONDO 2 Charging System. Includes adapter for US, UK, EU, and AU.

32734



Speech Processor Test Device

Quickly check if RONDO 2 is functioning by placing the flat side of the processor over the device. LED indicator lights up when RONDO 2 microphone is detecting sounds and coil is transmitting signal.

04048



Audio Processor Travel Case

Convenient zippered carrying case for audio processor and accessories.

33722



All-Purpose Backpack

High-quality backpack for carrying accessories and other items.

35280




Daily Case (Medium)

Portable case for carrying RONDO 2 and accessories.

34706

Stay Connected

Whether phones, music players, or tablets—RONDO 2 connects seamlessly to handheld devices and Bluetooth neckloops to provide optimal hearing, even in challenging environments.





Mini Battery Pack for RONDO 2

The Mini Battery Pack powers RONDO 2 for up to 35 hours on a single AAA battery. It can also run for up to 16 hours on a single rechargeable DaCapo PowerPack. The Mini Battery Pack is necessary to connect wireless receivers, FM systems, or direct audio input cables to RONDO 2.



RONDO 2 Mini Battery Pack Kit

Mini Battery Pack and Cable for RONDO 2, plus accessories.	36205
1 Mini Battery Pack & RONDO 2 Cable	35317
1 AAA Battery	08083
1 Mini Battery Pack Accessories	32579

Mini Battery Pack & RONDO 2 Cable

Mini Battery Pack and Cable for RONDO 2.	35317
1 RONDO 2 Mini Battery Pack Cable with Clip	31117
1 Mini Battery Pack	05880

Neckloops

Telecoil neckloops are a simple and convenient way to stream audio wirelessly to RONDO 2. You can connect to your phone through Bluetooth or plug into your favourite hand-held devices.



Artone 3 MAX Neckloop

Integrated Bluetooth and dual microphones enable hands-free phone calls and music streaming.

31776

Geemarc iLoop+ Neckloop

Connects to any standard audio headphone jack with a 3.5 mm plug.

08320

Direct Audio Input Cables

Connect RONDO 2 to audio devices, like phones, music players, and tablets. The 3.5 mm jack at one end of the cable plugs into any standard headphone port. The 3-pin jack at the other end connects to the Mini Battery Pack.



Unilateral Mixed Source

Provides a 50/50 mix of sound from the audio device and the RONDO 2 microphone, depending on settings. For use with one audio processor.

04438



Bilateral Mixed Source

Provides a 50/50 mix of sound from the audio device and the RONDO 2 microphone, depending on settings. For use with two audio processors.

04440



Unilateral External Source

Provides a 90/10 mix of sound, with more sound from the audio device, depending on settings. For use with one audio processor.

04439



Bilateral External Source

Provides a 90/10 mix of sound, with more sound from the audio device, depending on settings. For use with two audio processors.

04441

Additional Fixation

These easy-to-use accessories help keep RONDO 2 right where it should be. They provide additional security during sports and other activities.



Sports Headband

Soft and breathable headband with specially-designed pocket for RONDO 2—to keep it in place and protect it from perspiration and moisture during sports.

XS (<51 cm)	51403
S (51-55 cm)	51404
M (56-59 cm)	51405
L (>59 cm)	51406



Clothes & Hair Attachment Clip Set

Set containing both attachment clips	35318
Clothes Attachment Clip	35523
Hair Attachment Clip	35522



Fixation Loop for Glasses

Soft loop to attach RONDO 2 to the arm of glasses.	32239
--	-------



RONDO 2 Leash

8 cm Leash (Single)	33517
8 cm Leash (4 Pack)	36531
25 cm Leash (Single)	33518
25 cm Leash (4 Pack)	36532



GripWear

Reusable adhesive pad which is applied to the bottom of RONDO 2 for better fixation. Includes 12 pads per package.	33686
--	-------



Activity Clip

Small magnetic device to help keep RONDO 2 securely in place. Please note: not compatible with the SYNCHRONY implant.	51369
---	-------



WaterWear for RONDO 2

WaterWear is a fully waterproof accessory cover for RONDO 2 that can be used in freshwater or saltwater, as well as chlorinated water. Each package contains 3 reusable covers and 9 single-use adhesive strips for a total of 9 uses.



WaterWear for RONDO 2

Includes 3 covers and 9 adhesive strips for 9 total uses.

34747



RONDO 2 Clinical Kits

RONDO 2 Clinical Kits are designed for clinicians who are fitting RONDO 2. There's kits for programming RONDO 2, as well as kits for testing the function of the microphone.



RONDO 2 Clinical Kits

RONDO 2 Clinical Kit

The RONDO 2 Clinical Kit includes a RONDO 2 and charging pad for your clinic. This allows you to easily fit recipients while simultaneously charging their RONDO 2. It also includes the RONDO 2 programming cable and other accessories. 36200

- | | |
|-------------------------------------|----------------------------------|
| 1 RONDO 2 Processor Unit SP | 1 MAESTRO 7.0 Activation Key |
| 1 Charging Pad with USB Cable | 1 USB Power Supply |
| 1 Magnet Exchange Tool SUP | Supporting Manuals and Documents |
| 1 MAX Prog. Cable RONDO 2 with Clip | |

RONDO 2 Basic Clinical Kit

The RONDO 2 Basic Clinical Kit includes a charging pad, RONDO 2 programming cable, and other accessories. It does not include a RONDO 2 audio processor. 36201

- | | |
|-------------------------------------|----------------------------------|
| 1 Charging Pad with USB Cable | 1 USB Power Supply |
| 1 Magnet Exchange Tool SUP | Supporting Manuals and Documents |
| 1 MAX Prog. Cable RONDO 2 with Clip | |
| 1 MAESTRO 7.0 Activation Key | |

RONDO 2 MTD Kit

The RONDO 2 MTD Kit is designed to let you test the microphone function of RONDO 2. It includes the microphone test device, the Mini Battery Pack & Cable for RONDO 2, and other accessories. 36202

- | | |
|----------------------------|----------------------------------|
| 1 Microphone Test Device | 1 MBP MTD Connecting Cable |
| 1 MTD Headphones | 1 MBP & RONDO 2 Cable |
| 1 AAA Battery, for Mini BP | Supporting Manuals and Documents |
| 1 Battery 1.5V R6 for MTD | |

RONDO 2 Basic MTD Kit

The RONDO 2 Basic MTD Kit is designed for clinicians who already have a microphone test device. The Basic MTD Kit includes the Mini Battery Pack & Cable for RONDO 2, and other accessories. 36203

- | |
|----------------------------------|
| 1 AAA Battery for Mini BP |
| 1 MBP MTD Connecting Cable |
| 1 MBP & RONDO 2 Cable |
| Supporting Manuals and Documents |

RONDO 2 Magnet Kit

Kit for clinicians with additional RONDO 2 magnets of different strengths. Includes Magnet Exchange Tool SUP for changing RONDO 2 magnets. 36204

- | | |
|--------------------|----------------------------|
| 2 RONDO 2 Magnet 1 | 2 RONDO 2 Magnet 1S |
| 3 RONDO 2 Magnet 2 | 4 RONDO 2 Magnet 2S |
| 2 RONDO 2 Magnet 4 | 2 RONDO 2 Magnet 4S |
| | 1 Magnet Exchange Tool SUP |





Get In Touch

Where to Order

To order parts and accessories for RONDO 2,
contact your local MED-EL representative.

Let us know if you have any questions—
we look forward to hearing from you!



medel.com

MED-EL Offices Worldwide

AMERICAS

Argentina

medel@medel.com.ar

Canada

officecanada@medel.com

Colombia

office-colombia@medel.com

Mexico

office-mexico@medel.com

United States

implants@medelus.com

ASIA PACIFIC

Australia

office@medel.com.au

China

office@medel.net.cn

Hong Kong

office@hk.medel.com

India

implants@medel.in

Indonesia

office@id.medel.com

Japan

office-japan@medel.com

Malaysia

office@my.medel.com

Philippines

office@ph.medel.com

Singapore

office@sg.medel.com

South Korea

office@kr.medel.com

Thailand

office@th.medel.com

Vietnam

office@vn.medel.com

EMEA

Austria

office@at.medel.com

Belgium

office@be.medel.com

Denmark

nordic@medel.com

Finland

office@fi.medel.com

France

office@fr.medel.com

Germany

office@medel.de

Italy

ufficio.italia@medel.com

Netherlands

office@nl.medel.com

Norway

norge@medel.com

Portugal

office@pt.medel.com

South Africa

customerserviceZA@medel.com

Spain

office@es.medel.com

Sweden

nordic@medel.com

Switzerland

office-ch@medel.com

United Arab Emirates

office@ae.medel.com

United Kingdom

customerservice@medel.co.uk

26544 r1.0