

RONDO[®] Audio Processor

Product Catalogue



hearLIFE

RONDO

Audio Processor

The RONDO Audio Processor has a revolutionary single-unit design, combining the control unit, battery pack and coil into one compact device. RONDO is worn directly over the implant, leaving nothing on the ear. This makes it comfortable to wear, particularly for those with glasses, and it can be discreetly worn under the hair. And to give every MED-EL user these benefits, RONDO is compatible with all currently available MED-EL cochlear implants.



RONDO Base Kit



RONDO Base Kit

The RONDO Base Kit includes everything necessary to get started with RONDO. Please note, packages are subject to change and content may vary. Additional products, such as the Mini Battery Pack Kit, need to be ordered separately.

- | | |
|--|---------------------------------------|
| 1 RONDO Control Unit | 1 Speech Processor Test Device |
| 1 RONDO Battery Pack | 2 Attachment Clips (Clothes and Hair) |
| 2 RONDO Magnets
(2 and 3 / 2S and 3S) | 2 Fixation Loops for Glasses |
| 4 RONDO Microphone Covers | 1 GripWear |
| 2 Protectors (Silicone Cover) | 1 Carrying Bag |
| 1 FineTuner With Battery | 1 Travel Case |
| 2 Zinc-Air Battery Packages | 1 USB Stick |
| 1 Electrical Drying System | Supporting Manuals and Documents |

Anthracite	08798	Ebony	08800
Creme	08801	Nordic Grey	08799

RONDO Main Components



Control Unit

Includes Microphone Cover. The audio processor will need to be refitted if this component is replaced.

	Anthracite	08284		Ebony	08477
	Creme	08475		Nordic Grey	08476

Battery Pack



For use with 3 zinc-air size 675 batteries.

	Anthracite	08286		Ebony	08480
	Creme	08478		Nordic Grey	08479



Magnet SYNCHRONY

For SYNCHRONY implant only.

	1S—Soft	31005
	2S—Standard	31006
	3S—Strong	31007
	4S—Super Strong	31008
	5S—Max Strength	34365





Magnet

For all implants except SYNCHRONY.

	1—Soft	08280
	2—Standard	08281
	3—Strong	08282
	4—Super Strong	08283
	5—Max Strength	34367

Microphone Cover

Replacement microphone covers can be ordered individually or in sets of 4.

	1 Cover	4 Covers
	08024	32602
	08505	32603
	08518	32605
	08531	32604

RONDO Additional Components



FineTuner Remote Control

The FineTuner remote is used to change settings such as volume, program, and connectivity settings. FineTuner is compatible with all current MED-EL audio processors. One FineTuner can be used to control both processors for bilateral patients. Includes 1 CR2025 battery.

03132



Dry-Star UV

Daily drying helps keep RONDO in optimal working condition, and UV light sanitises audio processor. Powered by included AC adapter.

32338

FineTuner Battery

CR2025

03142



Speech Processor Test Device

Quickly check if RONDO is functioning by placing the flat side of the processor over the device. LED indicator lights up when RONDO microphone is detecting sounds and coil is transmitting signal.

04048



RONDO Protector

Silicone cover to protect RONDO from fall damage and other minor impacts.

08014



Audio Processor Travel Case

Convenient zippered carrying case for audio processor and accessories.

33722



Audio Processor Batteries

Package of 6 mercury-free, high-power, zinc-air size 675 batteries.

09564



Carrying Bag


Reusable MED-EL bag for carrying audio processor kit and other accessories.

32746



Stay Connected

Whether phones, music players, or tablets—RONDO connects seamlessly to handheld devices and Bluetooth neckloops to provide optimal hearing, even in challenging environments.





Mini Battery Pack for RONDO

The Mini Battery Pack powers RONDO for up to 37 hours on a single AAA battery. It can also run for up to 16 hours on a single rechargeable DaCapo PowerPack.



RONDO Mini Battery Pack Kit, 27 cm

1 RONDO Mini Battery Pack Cable (27 cm)	see below
1 Mini Battery Pack (Anthracite)	05880
1 AAA Battery	08083
1 Pouch	08246
1 Pouch Clip	08247
1 Pouch Lanyard	08248
1 User Manual (English)	06302

Anthracite Kit Cable	09983 08274
Creme Kit Cable	09989 08481

Ebony Kit Cable	09987 08483
Nordic Grey Kit Cable	09985 08482

RONDO Mini Battery Pack Kit, 100 cm

1 RONDO Mini Battery Pack Cable (100 cm)	see below
1 Mini Battery Pack (Anthracite)	05880
1 AAA Battery	08083
1 Pouch	08246
1 Pouch Clip	08247
1 Pouch Lanyard	08248
1 User Manual (English)	06302

Anthracite Kit Cable	09984 08275
Creme Kit Cable	09990 08484

Ebony Kit Cable	09988 08486
Nordic Grey Kit Cable	09986 08485

Neckloops

Telecoil neckloops are a simple and convenient way to stream audio wirelessly to RONDO. You can connect to your phone through Bluetooth or plug into your favourite hand-held devices.



Artone 3 MAX Neckloop

Integrated Bluetooth and dual microphones enable hands-free phone calls and music streaming.

31776



Geemarc iLoop+ Neckloop

Connects to any standard audio headphone jack with a 3.5 mm plug.

08320

Direct Audio Input Cables

Connect RONDO to audio devices, like phones, music players, and tablets. The 3.5 mm jack at one end of the cable plugs into any standard headphone port. The 3-pin jack at the other end connects to the RONDO Mini Battery Pack.



Unilateral Mixed Source

Provides a 50/50 mix of sound from the audio device and the RONDO microphone, depending on settings. For use with one audio processor.

04438



Bilateral Mixed Source

Provides a 50/50 mix of sound from the audio device and the RONDO microphone, depending on settings. For use with two audio processors.

04440



Unilateral External Source

Provides a 90/10 mix of sound, with more sound from the audio device, depending on settings. For use with one audio processor.

04439



Bilateral External Source

Provides a 90/10 mix of sound, with more sound from the audio device, depending on settings. For use with two audio processors.

04441

Additional Fixation

These easy-to-use accessories help keep RONDO right where it should be. They provide additional security during sports and other activities.



Sports Headband

Soft and breathable headband with specially-designed pocket for RONDO—to keep it in place and protect it from perspiration and moisture during sports.

XS (<51 cm)	51403
S (51-55 cm)	51404
M (56-59 cm)	51405
L (>59 cm)	51406



Clothes Attachment Clip

Small clip that secures RONDO discreetly to clothing.

08795



Hair Attachment Clip

Small clip that secures RONDO discreetly to hair.

30447

Clothes & Hair Attachment Clip Set

Set containing both attachment clips

32879



GripWear

Reusable adhesive pad which is applied to the bottom of RONDO for better fixation. Includes 12 pads per package.

33686



Activity Clip

Small magnetic device to help keep RONDO securely in place. Please note: not compatible with the SYNCHRONY implant.

51369



Fixation Loop for Glasses

Soft loop to attach RONDO to the arm of glasses.

32239



WaterWear for RONDO

WaterWear is a fully waterproof accessory cover for RONDO that can be used in freshwater or saltwater, as well as chlorinated water. Each package contains 3 reusable covers and 9 single-use adhesive strips for a total of 9 uses. Each WaterWear cover gives waterproof protection for up to 16 hours.

Only alkaline, silver-oxide or rechargeable 675 batteries can be used with WaterWear.



WaterWear for RONDO

Includes 3 covers and 9 adhesive strips for 9 total uses.

30240



Alkaline Batteries LR44

Compatible with WaterWear, provides up to 8 hours of battery life.

32328



Get In Touch

Where to Order

To order parts and accessories for RONDO, contact your local MED-EL representative.

Let us know if you have any questions—we look forward to hearing from you!



medel.com

MED-EL Offices Worldwide

AMERICAS

Argentina

medel@medel.com.ar

Canada

officecanada@medel.com

Colombia

office-colombia@medel.com

Mexico

office-mexico@medel.com

United States

implants@medelus.com

ASIA PACIFIC

Australia

office@medel.com.au

China

office@medel.net.cn

Hong Kong

office@hk.medel.com

India

implants@medel.in

Indonesia

office@id.medel.com

Japan

office-japan@medel.com

Malaysia

office@my.medel.com

Philippines

office@ph.medel.com

Singapore

office@sg.medel.com

South Korea

office@kr.medel.com

Thailand

office@th.medel.com

Vietnam

office@vn.medel.com

EMEA

Austria

office@at.medel.com

Belgium

office@be.medel.com

Denmark

nordic@medel.com

Finland

office@fi.medel.com

France

office@fr.medel.com

Germany

office@medel.de

Italy

ufficio.italia@medel.com

Netherlands

office@nl.medel.com

Norway

norge@medel.com

Portugal

office@pt.medel.com

South Africa

customerserviceZA@medel.com

Spain

office@es.medel.com

Sweden

nordic@medel.com

United Arab Emirates

office@ae.medel.com

United Kingdom

customerservice@medel.co.uk

23847 r3.0